


Chiropractic

DIGITAL RADIOGRAPHIC FLOOR/CEILING
MOUNTED SYSTEM PACKAGE

Universal

VIEWWORKS



Package Overview

Universal Chiropractic Systems uniquely combine economy and flexibility in a system designed to meet specialized and individual professional radiographic needs. This floor/ceiling system is thoroughly at home in private, group, or clinical applications. It provides full radiographic capabilities, from skull to weight bearing studies.

Standard Configuration

Vieworks DR System

- (1) 17x17 Gadolinium Oxysulfide
Wired cassette-sized detector

Generator

- AN-Series 30kW
with DR Integration

VS100 Vertical Stand

- Manual Cassette Tray
- 103 lp 10:1 grid

X-ray Tube

- Toshiba E7239

Collimator

- Ralco Manual

** Additional options and upgrades available*