

# DELWORKS **FIT**

## PORTABLE WORKSTATION

The DELWORKS FIT portable tablet workstation aims to maximize the portability and efficiency of the DELWORKS software - utilizing it anytime, anywhere.

All of the excellent user-friendly features of the DELWORKS Software are beautifully translated to an able-bodied tablet PC, capable of withstanding even the harshest radiographic environments. This handheld workstation features wireless connectivity, plenty of storage space, and a powerful Intel® Core™ i5 processor to perform rapid and productive patient side studies, without limitations.

## Key Highlights

### Full HD Resolution

-12.5" wide viewing angle display

### Long-lasting Battery

-Over 9 hours of battery life

### Dual Network Interface Cards (NIC)

-For increased range and stability with Hospital network and dedicated detector connectivity

### Wirelessly Transfer Images Directly to PACS

### Ultra-light and Portable

- Weighs only 2.9 lbs. (1.34 kg)

### Dockable Unit

-Allowing charging, transferring, and screen mirroring

### Rugged, Protective Design

-Tested to handle multiple 4' drops, water, dust, and splash resistance




# Technical Specifications

\* All specifications subject to change without notice

## DELWORKS Fit Portable Workstation

Display:	Full-HD resolution (1920 x 1080) / 12.5" wide viewing angle screen
Processor:	Intel® Core™ i5 (2.3GHz)
Storage:	250GB solid state drive (SSD)
Memory (RAM):	8GB
Operating system:	Windows 10 Pro 64-bit
Keyboard/Input:	Capacitive 10-point touchscreen
Battery:	Battery life over 9 hours / approximately 2.5 hours charge time
Weight:	2.9 lbs. (134 kg)
Dimensions:	12.93" x 8.17" x .75"
Certifications:	MIL-STD-810G tested (multiple 4' drops) / IP54 (water, dust, and splash resistant)
Wifi:	Dual-band wireless-ac 8260 wi-fi plus bluetooth® 4.2 Wireless-n 150Mbps nano usb adapter
Dockable:	Yes



[www.delmedical.com](http://www.delmedical.com)

**New York**  
28 Calvert Street  
Harrison, NY 10528  
(Tel) 800.261.9808  
(Fax) 914.835.6111

**Illinois**  
241 Covington Drive  
Bloomington, IL 60108  
(Tel) 800.800.6006  
(Fax) 847.288.7011

**DEL MEDICAL**

#9000-DW-FITWorkstation 7.19